



Titin Antibody

Titin Antibody (TitinAb) ELISA Test Kit[†]

KRONUS[®], a leading provider of specialized autoimmune diagnostic test kits, is pleased to announce the availability of a new ELISA test kit for the measurement of antibodies to titin (TitinAb). Titin is a protein of the striated muscles with an extremely high molecular weight. The immunogenic region of titin is located on a 30kD protein fragment.

Antibodies against this fragment presumably cross-react with the epitopes of the acetylcholine receptors making their measurement potentially useful in select indications.

This new ELISA utilizes recombinant MGT30 peptide providing for the accurate and specific determination of titin antibodies.

To obtain additional information on **KRONUS'** unique and progressive product line, please call us toll-free at **800 4 KRONUS** or visit us at our web site at www.kronus.com.



Your Source for Sensitive Autoimmune Diagnostics

[†] For Research Use Only. Not for use in diagnostic procedures.

ISO 13485 : 2003 QMS Certified
201001